








*If you need further assistance, please call
your local Custom Medical Solutions
Representative or Customer Service at
1-866-350-5640*

Features



TROUBLESHOOTING ALARMS

ALARM TYPE	CONDITION	ACTION/RESOLUTION	USER TIPS
<p style="text-align: center;">▶ TARGET TIME OUT</p> <p style="text-align: center;">(LEAKAGE)</p>	<p>There may be a leak in the dressing or the canister</p> <p>This alarm will sound when selected therapy pressures are not being reached or maintained.</p>	<ol style="list-style-type: none"> 1. Ensure dressing/canister lure lock connector is properly locked 2. Ensure canister is fully engaged 3. Identify air leaks in the dressing 	<ul style="list-style-type: none"> • Identify air leaks with stethoscope • Use additional drape to secure air leaks in dressing • Make sure dressing is fully collapsed • Make sure therapy unit and tubing is level with dressing tubing
<p style="text-align: center;">▶ BLOCKAGE</p> <p><i>WARNING: According to clinician instructions, replace XLR8 dressing with wet to dry or alternate dressing if therapy is interrupted or off for more than two hours</i></p>	<p>Unit has determined a blockage between the wound dressing and the canister</p>	<ol style="list-style-type: none"> 1. Check tubing for closed clamps, kinks, crimps, or blockages. 2. Inspect to ensure a quarter-sized (2.5cm) hole has been cut in the drape. 	<ul style="list-style-type: none"> • A smaller hole or a slit may cause blockage, pressure fluctuations, and other functionality issues. • Ensure port pad is located on a flat area of the body, avoiding skin folds
<p style="text-align: center;">▶ LOW BATTERY</p> <p style="text-align: center; color: red;">CRITICALLY LOW →→</p>	<p>Battery has 2% to 25% capacity</p>  <p>Battery has 0% to 2% capacity or 20 minutes of operation</p> 	<ol style="list-style-type: none"> 1. Recharge by plugging the XLR8 power adapter into the wall outlet 2. Ensure that the power cord is securely placed in the XLR8 charge port. 3. Lightning bolt icon indicates that the unit is receiving power from the adapter and is charging from the wall receptacle... 	<ul style="list-style-type: none"> • May take up to 3 hours to fully recharge • Unless the patient is ambulating or transferring we recommend that the power cord stays plugged into the power adapter 
<p style="text-align: center;">▶ CANISTER FULL</p> <p><i>WARNING: According to clinician instructions, replace XLR8 dressing with wet to dry or alternate dressing if therapy is interrupted or off for more than two hours.</i></p>	<p>Unit has detected that canister is full and must be replaced</p>	<p>Change canister and restart therapy</p>	<ul style="list-style-type: none"> • If minimal exudate is noted in canister the filter may have gotten wet. Canister will need to be changed to reset system • Keep therapy unit upright at all times to avoid unnecessary changes <p>(Correct) (Incorrect)</p>  

For additional support call Custom Medical Solutions: 1-866-350-5640